

3D Printed Screw & Rod System

VENTOTENE


Straight and
Pre-Bended Rods




STRAIGHT AND PRE-BENDED RODS

52-20/21/22/24 YY.S

52-60YY.S



Code	Diameter xx (mm)	Length yy (mm)	Code	Diameter xx (mm)	Length yy (mm)
52-1040.S	5	40	52-2040.S	5,5	40
52-1050.S	5	50	52-2050.S	5,5	50
52-1060.S	5	60	52-2060.S	5,5	60
52-1070.S	5	70	52-2070.S	5,5	70
52-1080.S	5	80	52-2080.S	5,5	80
52-1090.S	5	90	52-2090.S	5,5	90
52-1100.S	5	100	52-2100.S	5,5	100
52-1110.S	5	110	52-2110.S	5,5	110
52-1120.S	5	120	52-2120.S	5,5	120
52-1140.S	5	140	52-2140.S	5,5	140
52-1160.S	5	160	52-2160.S	5,5	160
52-1180.S	5	180	52-2180.S	5,5	180
52-1200.S	5	200	52-2200.S	5,5	200
52-1450.S	5	450	52-2450.S	5,5	450



Code	Diameter xx (mm)	Length yy (mm)	Code	Diameter xx (mm)	Length yy (mm)
52-5040.S	5	40	52-6035.S	5,5	35
52-5050.S	5	50	52-6040.S	5,5	40
52-5060.S	5	60	52-6045.S	5,5	45
52-5070.S	5	70	52-6050.S	5,5	50
52-5080.S	5	80	52-6055.S	5,5	55
52-5090.S	5	90	52-6060.S	5,5	60
52-5100.S	5	100	52-6065.S	5,5	65
52-5110.S	5	110	52-6070.S	5,5	70
52-5120.S	5	120	52-6075.S	5,5	75
52-5140.S	5	140	52-6080.S	5,5	80
52-5160.S	5	160	52-6090.S	5,5	90
52-5180.S	5	180			

Materials

- Titanium Ti6Al4V ELI

**Dimensional tolerance $\pm 0,1$ mm

Tsunami Medical S.r.l.
Via XXV Aprile n. 22
41037 Mirandola (MO) Italy
Tel. +39.0535.38397
Fax. +39.0535.38399
www.tsunamimedical.com
info@tsunamimed.com

CE 0051
UNI CEI EN ISO 13485

Made in Italy

BR_RODS_REV0_20240513